



SF40A

Impact Printer



OHAUS Value in a Full-Featured, Portable Printer

The SF40A combines high performance and multiple functions to offer OHAUS quality and durability in a portable printer. The SF40A offers advanced functionality in an easy-to-use and set up printer at an exceptional value. This versatile printer is ideal for use in conjunction with OHAUS balances and scales in laboratory, industrial, and educational applications to ensure accurate documentation of weighing results.

Unique Features:

- The statistics and totalization modes and real time clock complement OHAUS weighing instruments' application modes and offer functionality that assists the reporting processes.
- SF40A helps support the efficiency of the weighing and recording processes. RS232 interface and PC configuration software make setting up and using the printer fast and simple.
- With the ability to print in 13 different languages, SF40A can print in the language of the connected device to support universal operation.

SF40A Impact Printer

Model	SF40A
Number of Columns	40 columns
Print Speed	Approximately 1.2 lines per second
Application Modes	Totalization, Statistics (statistics printout includes number of samples, average value, standard deviation, minimum, maximum, min-max differential, and total sample weight sum)
Paper Roll (W x Dia.)	Standard paper 57.5 mm x Ø 60 mm max.
Power Requirements	AC adapter Input: 100 – 240 VAC, 50 - 60 Hz AC adapter output: 12 VDC, 2.5 A
Interface	RS232C
Ambient Temperature	0 °C – 40 °C
Printer Dimensions (W x L x H)	140 x 205 x 75 mm
Net Weight	740 g (Paper roll included)
Shipping Weight	0.8 kg

Other Standard Features

Status indicator, 40 characters per line, perform zero and tare operations directly from the printer. Full functionality with Adventurer™, Explorer™, Navigator™, MB Series Moisture Analyzers and 'T' Series indicators. Statistics and Totalization not supported on CKW Bench Scale and Scout™ Series. Multi-Language Printing: English, Spanish, German, Italian, French, Chinese, Korean, Japanese, Russian, Polish, Portuguese, Hungarian and Czech.


Accessories

Paper (Included).....	12120799
Paper Rotary Axis (Included)	30063920
Ink Ribbon (Included)	30529322
Adhesive Paper Rolls (3 pcs)	11600388
Printer Cable RS9 M/F (Included)	30529323
Cable Adapter (Included).....	30059316
Printer Cable RS9 F/Bare Wire.....	80500552
Power Adapter, SF40A (Not Included)	30529321

Compliance

- **Product Safety:** IEC/EN 61010-1
- **Electromagnetic Compatibility:** IEC/EN 61326-1 Class B, Basic environments; FCC Part 15 Class A; Canada ICES-003 Class A
- **Compliance Marks:** CE, RCM, UKCA, CSA
- **Other:** WEEE, RoHS

Outline Dimensions



OHAUS Europe GmbH

Heuwinkelstrasse 3,
8606 Nänikon,
Switzerland

e-mail: ssc@ohaus.com
Tel: 0041 22 567 53 19
e-mail: tsc@ohaus.com
Tel: 0041 22 567 53 20

www.ohaus.com

The management system
governing the manufacture of
this product is ISO 9001:2015
certified.

